MANUAL DE INSTRUCCIONES

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PISCIS

Manual de instrucciones **ES**

Instruction manual **EN**

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Manuale D'Istruzioni **IT**

Anweisungen **DE**

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Bruksanvisning **SWE**

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NOTICE

Notice to the user and/or patient: any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is established.

1. INSTRUCTIONS

1. The device should only be used by one user.
2. Always use a slower speed when performing the following actions: ascending, descending or traversing a slope or inclined terrain, as well as on uneven terrain, ramps and soft surfaces.
3. Keep away from flammable objects while the device is charging.

2. CONTRAINDICATIONS



1. Never attempt to ascend or descend a slope on slippery surfaces or where there is a danger of skidding.
2. Do not drive on a slope that exceeds 6 degrees.
3. Avoid shifting the centre of gravity when turning left or right.
4. Direct exposure to rain or moisture will cause the scooter to malfunction and may cause the scooter to rust.
5. Never put the scooter in neutral when on slopes.
6. Do not attempt to modify or tamper with the mechanism or structure of the scooter.
7. Do not use the scooter on stairs or escalators.
8. Do not use the device while it is charging.

3. GETTING STARTED

1. To start using the scooter, carry out the following steps:
 - Remove the rear casing and connect the battery.
 - Check that the immobiliser lever is in D mode.
 - Check that the speed control is in low speed (turtle) mode.
2. Stand on the scooter and start the scooter with the key. The battery indicator will light up. Confirm that the self-diagnostic indicator is not flashing.
3. By placing your hands on the handlebars, you should comfortably reach the control levers, which are operated with your thumbs. The right-hand lever activates the scooter to move forwards and the left-hand lever moves the scooter backwards. If you release both levers, the scooter will stop.
4. Practice driving where there are no obstacles, starting at a slow speed.
5. If the scooter stops, locate the circuit breaker (A). Press it and test the scooter again.

4. IDENTIFICATION OF PARTS

Before using the scooter, it is important that you are familiar with the parts of the scooter and the operation of the controls and controls.

5. FUNCTIONS OF THE PARTS

Control panel

- 1. Speed control:** The tortoise icon indicates a slow speed level. The hare icon indicates a fast speed level. Before starting to use the scooter, select the slow speed level (turn the wheel towards the tortoise icon). Increase the speed slightly when driving uphill and reduce the speed when driving downhill.
- 2. Self-diagnostic warning indicator:** A flashing light indicates that there is a problem with the scooter.
- 3. Battery level indicator:** When all five stripes are lit, the battery is fully charged. As they fade, the battery is running low.
- 4. Throttle:** Pressing the right throttle will move the scooter forward. Pressing the left accelerator moves

the scooter backwards. Releasing them activates the automatic brake. Depending on how hard you press the accelerator, it will move at a faster or slower speed. Keep a safe distance to avoid any danger.

Handlebar base

1. Adjusting the steering handlebars: Push the lever down to adjust. Pull to lock the tiller angle adjuster at a comfortable angle.

Under the seat

1. Seat adjustment lever: Pull the lever backwards to rotate the seat; push it forward to lock the seat position.

2. Armrest adjustment: Loosen the two thumbwheels to adjust the armrest; tighten them again to lock in the desired position.

Rear part

1. Rear cover: Used to protect the battery from damage.

2. Anti-tip wheels: Used to prevent the scooter from tipping over.

3. Starter lever: The starter lever can be used to switch the scooter into manual mode (N) or automatic mode (D). Manual mode allows the user to move the scooter without using the battery, i.e. by pushing the scooter, while automatic mode allows the user to both drive the scooter and keep it stationary when parked.

Indicators on the lever will help you determine which mode the scooter is in. Whenever you want to change from manual to automatic mode, you must first switch the scooter off. You must then change mode using the lever. Then switch on the scooter to start using it.

Observe the following instructions:

- Release the throttle and allow the scooter to come to a complete stop before shifting from forward to reverse and vice versa.
- Use the scooter only in areas where it is safe to walk.
- Use the low speed mode for downhill slopes, reversing, ramps, kerbs or uneven surfaces. For ascending slopes, you may need to use a higher speed.
- For descending steep slopes, set the scooter to low speed mode.

6. ELECTROMAGNETIC INTERFERENCE WARNINGS

The operation of power wheelchairs and motorized scooters may be affected by electromagnetic interference (EMI). This is the result of interference from electromagnetic (EM) energy waves emitted by sources such as radio stations, television stations, amateur radio transmitters, two-way radios and mobile phones. Interference (from radio wave sources) could cause the braking system to release unintentionally, cause the scooter to move by itself or to move in an undesirable direction. They may also cause irreparable damage to the scooter's control system.

1. Do not operate portable transceivers (transceivers), such as CB (Citizens Band) radios. Do not activate personal communication devices, such as mobile phones, while the scooter is switched on.

2. Be aware of the existence and location of any nearby transmitters, such as radio or television stations,

and avoid getting too close to them.

3. In the event of unintended movement or brake release, switch the scooter off.
4. Please note that fitting accessories or components to the scooter or modifying the scooter in any way may make it more susceptible to the effect of electromagnetic interference.

The immunity level of 20 V/m is generally an affordable and practical level of protection against electromagnetic interference (as of May 1994). Remember that the higher the immunity level, the higher the level of protection. This product has an immunity level of at least 10 V/m.


7. BATTERY CHARGING

For the first time, it is advisable to let the battery charge overnight and it should be recharged after each day's use.

1. Insert the battery cable into the connector.
2. Plug the other end of the cable into an electrical outlet.
3. Keep the battery charged until the light turns green.
4. Once the light turns green, disconnect the battery from the connector.

If the green light flashes for more than 40 minutes, there is a problem with charging.

Check that the charger is correctly positioned and that the scooter is switched off. If the problem persists, reset the system by pressing the circuit breaker (A).

Charge the batteries after each use. Check the battery indicator before using the scooter again. If the scooter is not used for a period of time, the batteries may become discharged. Remember to charge the batteries at least once a month.

It is recommended that you charge the batteries fully each time so as not to damage the life of the batteries. Failure to do so may result in malfunctions.

8. DISASSEMBLING THE SCOOTER

The scooter fits easily into the boot or back seat of most cars for transportation. However, if you need to disassemble it, no tools are required. Just follow the reverse assembly instructions.


Here are a few brief guidelines to help you get your bearings:

1. Release the rotation lever to remove the seat.

2. Remove the rear cover of the scooter.
3. Remove the Velcro fastening from the battery.
4. Disconnect the power plugs from the battery.
5. Remove the two batteries.
6. Remove the basket and tilt the steering handlebar.
7. You should have 5 main parts after disassembling the scooter.

9. FUSE REPLACEMENT

1. Remove the front cover screws (R), the arrow indicates the direction (Caution, the screws have different sizes).
2. Next, open the fuse box (R1) and replace it with a back-up fuse (R2); note that the fuse must be of the same rating.

10. CLEANING AND MAINTAINING THE SCOOTER

Cleaning the scooter: If the scooter is dirty, use a damp or slightly soapy cloth to clean it. Do not wet it directly or rinse it, as it has electrical parts that should not come into direct contact with water. You can use car polish liquid.

All maintenance and repair work on the scooter must be carried out by a specialist.

The following areas require regular inspections:

- Tyre pressure should be between 30 and 50 psi.
- All electrical connections should be checked to ensure that they are securely connected.

TECHNICAL SPECIFICATIONS	
Total length	166 cm
Overall width (with armrests)	59 cm
Total height	95 cm
Front wheels	23 cm
Rear wheels	23 cm
Weight	67.5 kg
Maximum speed	8 km/h
Maximum supported weight	135 kg
Ground clearance	4 cm
Degree of inclination	12 degrees
Maximum kerb height	5 cm
Turning radius	141.5 cm
Brake	Electromechanical
Type of seat	Rotatable, lightweight and foldable
Seat width	46 cm
Engine size	450W, 3700 rpm
Batteries	(2) 12V, 36Ah
Battery pack weight	22.5 kg
Autonomy	34 km
Battery charger	5A external

11. PROBLEM DIAGNOSIS

Count the number of times the blue light on your scooter flashes to find out its fault.

NUMBER OF FLASHES	DESCRIPTION	MEANING
1	Low battery	The batteries are exhausted. - Recharge the battery.
2	Battery depleted	The batteries are exhausted. - Recharge the battery.

3	Battery overload fault	Battery voltage is too high. If the charger is connected, disconnect it. May have been caused by overcharging or overcharging or driving downhill. - When driving downhill, reduce speed. - Check its condition and connections.
4	Engine runtime failure	The engine running time has been exceeded. Switch the scooter off and leave it for a few minutes before switching it on again. - Check the motor, connections and wiring for short circuits. there is a short circuit. - Switch off the controller and allow a few minutes before switching back on.
5	Immobilisation brake failure	Immobiliser brake is activated or faulty: - Make sure that the starter lever is in automatic mode (D). mode (D). - If the flashing light does not come on until the throttle is used, check that there is no short circuit in the brake. - If the light appears as soon as the scooter is switched on, check that there is no open circuit fault. - Check the brake, its connections and wiring.
6	Throttle failure	The throttle is not in neutral when the scooter is switched on. - Make sure the throttle is in neutral mode before starting your scooter. - The throttle may need to be re-adjusted. - Check the throttle connections.
7	Speed failure	Throttle or throttle connections are faulty. - Check wiring for short circuits or open connections. open connections. - Throttle may not have been installed correctly.
8	Engine failure	The motor or its wiring is faulty. - Check the motor, its connections and wiring.
9	Another mistake	The controller may have an internal fault. - Check connections and wiring.

ES - Peso máximo 135 kg
EN - Maximum weight 135 kg
IT - Peso massimo 135 kg
FR - Poids maximum 135 kg
PT - Peso máximo 135 kg
DE - Maximales Gewicht 135 kg
NL - Maximaal gewicht 135 kg
SWE - Högsta vikt 135 kg
PL - Maksymalna waga 135 kg
DK - Maksimal vægt 135 kg

CHIEN TI ENTERPRISE CO., LTD
No. 13, Lane 227, Fu Ying Rd., Hsin Chuang District, New Taipei City,
Taiwan
sales@chienti.com.tw.

MEDNET EC-REP GMBH
Borkstrasse 10, 48163, Muenster, Germany
ecrep@medneteuropa.com

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GRUPO R. QUERALTÓ S.A.